



PRODUCT INFORMATION

MARTRON INC.[®] P O BOX 309; MONROE, NC 28111 PH: (704) 289-1934; FAX: (704) 283-7400

EMERGENCY - **CHEMTREC 800-424-9300**

REF # - MFC-007526

MARTRON MTC TEMPORIZER OPERATING INFORMATION FOR TRIVALENT CHROMIUM PLATING SOLUTIONS

INTRODUCTION

Martron MTC TempORIZER was developed to provide TVC users with a means of reducing the cosmetic effects of small amounts of organic contaminants that produce darkness in mid current density areas.

Martron MTC TempORIZER is economical to use since only small additions are normally required.

APPLICATION

Martron MTC TempORIZER - Should only be used when approved by *MARTRON INC.*

USE

The amount of **Martron MTC TempORIZER** to use will depend on the degree of contamination as indicated by mid current density darkness and/or streaking. As a starting point, the addition of 0.1 - 0.2 mL/L of **Martron MTC TempORIZER** should suffice in most cases. If contamination symptoms persist, additional amounts in increments of (0.1 mL/L) should be made. Additions totaling greater than (1.0 mL/L) are not recommended.